

Accelx™ Implant System Motor



Ref #AX200



AX200HP (Optic)
Max torque value 70 Ncm

Ease of Use

Allows the user to quickly choose between various functions with the touch of a button

9 Memory Programs

Easily select between 9 pre-saved programs.

Advanced Irrigation Pump

The super-quiet pump provides for easy set-up of the irrigation tubes

Foot Pedal Control

The foot pedal provides rotation direction of the handpiece, coolant flow, program selection and speed, allowing the user to be hands-free

Improved Contra-Angle, Optic Handpiece

The reduced angle head with built-in spray nozzle provides deeper access into the mouth during surgery

Powerful Optic Micromotor

The micromotor is equipped with powerful optic function which helps procedural accuracy



Powerful Optic
Micromotor



Foot Pedal

Ref #AIM92820LI,R0